

Color Guide Options

XCL-1 (Single Layer Teeth)

AvaDent Extreme Linked Material	Vita Classical	Remark
BN00	A1	or Lighter
BN10	A1-A2	
BN20	A2-A3	
BN30	A3-A3.5	
BN35	A3.5	or Darker
YW10	B1	or Lighter
YW20	B2	or Lighter
GY10	C1	or Lighter
GY20	C2	or Lighter
RD20	D2	or Darker
BL00		Pure White
BL30	OM3	



Part Number: XCL-1SHADEGUIDE

Printed Base Shades Available

TrueDent Original

TrueDent Light

TrueDent Dark

XCL-2 (Double Layer Teeth)

AvaDent Extreme Linked Material	Vita Classical
A1	A1
A2	A2
A3	A3
A3.5	A3.5
B1	B1
B2	B2
C1	C1
C2	C2
D2	D2
OM1	OM1
OM3	OM3

NextDent MFH	Vita Classical	Remark
Bleach	B1	or Lighter
N1	A1-A2	

Milled Base Shades Available

AvaDent Original

Lucitone 199 Original

AvaDent Light

Lucitone 199 Light

AvaDent Dark

Lucitone 199 Medium Dark

Lucitone 199 Dark Pink