

Upper Anteriors *(Curve width from distal canines)*

Dentsply Portrait Moulds

11G (52.0 mm)	21G (56.0 mm)	25G (51.5 mm)	43D (45.5 mm)	62D (48.5 mm)
11H (54.5 mm)	21J (58.0 mm)	31F (49.5 mm)	43F (49.0 mm)	62G (52.5 mm)
12E (50.0 mm)	21X (50.0 mm)	32E (49.0 mm)	45F (51.5 mm)	65G (52.0 mm)
12G (53.0 mm)	22C (44.5 mm)	32G (52.5 mm)	45H (54.0 mm)	74E (49.0 mm)
13E (47.5 mm)	22E (49.0 mm)	41J (58.0 mm)	52C (44.0 mm)	74H (54.0 mm)
21C (45.5 mm)	22G (51.0 mm)	42D (48.0 mm)	54F (49.0 mm)	75E (49.5 mm)
21D (47.0 mm)	22H (54.0 mm)	42F (50.0 mm)	55D (47.5 mm)	75G (53.5 mm)
21E (51.0 mm)	24F (48.0 mm)	42G (52.0 mm)	55F (50.5 mm)	

Ivoclar BlueLine Moulds

A11 (45.5 mm)	A24B (47.0 mm)	A44 (45.5 mm)	P12 (54.5 mm)	P3 (47.5 mm)
A12 (49.0 mm)	A25 (51.0 mm)	A54 (50.5 mm)	P13 (49.0 mm)	P4 (49.0 mm)
A13 (47.5 mm)	A26 (51.0 mm)	A56 (51.5 mm)	P14 (45.5 mm)	P5 (49.5 mm)
A14 (45.5 mm)	A27 (52.5 mm)	A66 (50.0 mm)	P15 (45.5 mm)	P6 (49.5 mm)
A15 (53.0 mm)	A32 (49.0 mm)	A68 (53.5 mm)	P16 (51.0 mm)	P7 (50.0 mm)
A16 (54.5 mm)	A36 (50.5 mm)	A69 (53.5 mm)	P17 (48.0 mm)	P8 (52.0 mm)
A17 (56.0 mm)	A37 (52.0 mm)	P1 (47.0 mm)	P18 (58.0 mm)	P9 (52.5 mm)
A22 (47.0 mm)	A41 (46.0 mm)	P10 (53.0 mm)	P2 (47.0 mm)	
A24 (46.5 mm)	A42 (48.0 mm)	P11 (53.0 mm)		

Ivoclar IvoStar Moulds

1 (40.4 mm)	4 (46.8 mm)	32 (45.0 mm)	35 (47.8 mm)	43 (48.4 mm)
T2 (43.0 mm)	5 (47.7 mm)	33 (45.9 mm)	41 (42.2 mm)	44 (47.3 mm)
3 (44.7 mm)	31 (42.1 mm)	34 (45.4 mm)	42 (45.0 mm)	45 (48.7 mm)

Candulor NFC+ Moulds

550 (41.8 mm)	556 (47.6 mm)	664 (46.6 mm)	772 (48.4 mm)	880 (48.0 mm)
552 (43.1 mm)	660 (45.3 mm)	666 (47.9 mm)	774 (49.2 mm)	882 (48.3 mm)
554 (44.6 mm)	662 (45.5 mm)	770 (47.7 mm)	776 (49.3 mm)	884 (51.1 mm)

Dentsply Portrait IPN



Ivoclar BlueLine



Ivoclar IvoStar



Candulor NFC+



IF PRINTED, SCAN QR CODE. IF YOU ARE ON MOBILE, CLICK QR CODE FOR LINK.

Upper Anteriors *(Curve width from distal canines)*

Candulor TCR Moulds

203 (45.4 mm)	263 (46.9 mm)	304 (47.8 mm)	342 (42.2 mm)
204 (48.7 mm)	264 (47.8 mm)	332 (40.4 mm)	343 (45.9 mm)
262 (43.0 mm)	303 (45.1 mm)	333 (45.0 mm)	452 (44.8 mm)

Lower Anteriors *(Curve width from distal canines)*

Dentsply Portrait Moulds

C (35.0 mm)	G (38.0 mm)	K1 (43.5 mm)	O (40.0 mm)	U (42.0 mm)
D (33.5 mm)	H (36.5 mm)	L (37.0 mm)	P (39.5 mm)	V (41.5 mm)
E (38.0 mm)	I (43.0 mm)	M (39.5 mm)	R (40.5 mm)	W (42.0 mm)
F (36.5 mm)	J (40.5 mm)	N (38.0 mm)	S (40.0 mm)	

Ivoclar BlueLine Moulds

A2 (37.0 mm)	A4 (36.0 mm)	A6 (38.5 mm)	A8 (42.0 mm)
A3 (36.0 mm)	A5 (38.5 mm)	A7 (40.0 mm)	A9 (40.5 mm)

Ivoclar IvoStar Moulds

11 (36.0 mm)	13 (32.9 mm)	15 (35.2 mm)
12 (34.2 mm)	14 (33.1 mm)	16 (37.6 mm)

Candulor NFC+ Moulds

990 (35.0 mm)	994 (38.0 mm)	998 (40.0 mm)
992 (36.7 mm)	996 (37.7 mm)	

Candulor TCR Moulds

53 (36.1 mm)	81 (33.2 mm)	83 (37.6 mm)
80 (32.8 mm)	82 (35.3 mm)	

Candulor TCR Moulds



AvaDent Tooth Mould Selection Ruler



AvaDent Mould and Shade Charts



IF PRINTED, SCAN QR CODE. IF YOU ARE ON MOBILE, CLICK QR CODE FOR LINK.