

# UPPER ANTERIORS (Curve width from distal canines)

### **Dentsply Portrait Moulds**

11G (52.0 mm)	21G (56.0 mm)	25G (51.5 mm)	43D (45.5 mm)	62D (48.5 mm)
11H (54.5 mm)	21J (58.0 mm)	31F (49.5 mm)	43F (49.0 mm)	62G (52.5 mm)
12E (50.0 mm)	21X (50.0 mm)	32E (49.0 mm)	45F (51.5 mm)	65G (52.0 mm)
12G (53.0 mm)	22C (44.5 mm)	32G (52.5 mm)	45H (54.0 mm)	74E (49.0 mm)
13E (47.5 mm)	22E (49.0 mm)	41J (58.0 mm)	52C (44.0 mm)	74H (54.0 mm)
21C (45.5 mm)	22G (51.0 mm)	42D (48.0 mm)	54F (49.0 mm)	75E (49.5 mm)
21D (47.0 mm)	22H (54.0 mm)	42F (50.0 mm)	55D (47.5 mm)	75G (53.5 mm)
21E (51.0 mm)	24F (48.0 mm)	42G (52.0 mm)	55F (50.5 mm)	

#### **Ivoclar BlueLine Moulds**

A11 (45.5 mm)	A22 (47.0 mm)	A36 (50.5 mm)	A66 (50.0 mm)	P13	(49.0 mm)	P3	(47.5 mm)
A12 (49.0 mm)	A24 (46.5 mm)	A37 (52.0 mm)	A68 (53.5 mm)	P14	(45.5 mm)	P4	(49.0 mm)
A13 (47.5 mm)	A24B (47.0 mm)	A41 (46.0 mm)	A69 (53.5 mm)	P15	(45.5 mm)	P5	(49.5 mm)
A14 (45.5 mm)	A25 (51.0 mm)	A42 (48.0 mm)	P1 (47.0 mm)	P16	(51.0 mm)	P6	(49.5 mm)
A15 (53.0 mm)	A26 (51.0 mm)	A44 (45.5 mm)	P10 (53.0 mm)	P17	(48.0 mm)	P7	(50.0 mm)
A16 (54.5 mm)	A27 (52.5 mm)	A54 (50.5 mm)	P11 (53.0 mm)	P18	(58.0 mm)	P8	(52.0 mm)
A17 (56.0 mm)	A32 (49.0 mm)	A56 (51.5 mm)	P12 (54.5 mm)	P2	(47.0 mm)	P9	(52.5 mm)

### Ivoclar IvoStar Moulds Candulor NFC+ Moulds

1	(40.4 mm)	34 (45.4 mm)	550 (41.8 mm)	770 (47.7 mm)
2	(43.0 mm)	35 (47.8 mm)	552 (43.1 mm)	772 (48.4 mm)
3	(44.7 mm)	41 (42.2 mm)	554 (44.6 mm)	774 (49.2 mm)
4	(46.8 mm)	42 (45.0 mm)	556 (47.6 mm)	776 (49.3 mm)
5	(47.7 mm)	43 (48.4 mm)	660 (45.3 mm)	880 (48.0 mm)
31	(42.1 mm)	44 (47.3 mm)	662 (45.5 mm)	882 (48.3 mm)
32	(45.0 mm)	45 (48.7 mm)	664 (46.6 mm)	884 (51.1 mm)
33	(45.9 mm)		666 (47.9 mm)	

### **Candulor TCR Moulds**

203 (45.4 mm)	332 (40.4 mm)
204 (48.7 mm)	333 (45.0 mm)
262 (43.0 mm)	342 (42.2 mm)
263 (46.9 mm)	343 (45.9 mm)
264 (47.8 mm)	452 (44.8 mm)
303 (45.1 mm)	
304 (47.8 mm)	



## LOWER ANTERIORS (Curve of width from distal canines)

### **Dentsply Portrait Moulds**

C (35.0 mm)	J (40.5 mm)	R (40.5 mm)
D (33.5 mm)	K1 (43.5 mm)	S (40.0 mm)
E (38.0 mm)	L (37.0 mm)	U (42.0 mm)
F (36.5 mm)	M (39.5 mm)	V (41.5 mm)
G (38.0 mm)	N (38.0 mm)	W (42.0 mm)
H (36.5 mm)	O (40.0 mm)	
(43.0 mm)	P (39.5 mm)	

### **Ivoclar BlueLine Moulds**

A2 (37.0 mm)	A7	(40.0 mm)
A3 (36.0 mm)	A8	(42.0 mm)
A4 (36.0 mm)	A9	(40.5 mm)
A5 (38.5 mm)		
A6 (38.5 mm)		

#### **Ivoclar IvoStar Moulds**

16 (37.6 mm)

11 (36.0 mm)	990 (35.0 mm)
12 (34.2 mm)	992 (36.7 mm)
13 (32.9 mm)	994 (38.0 mm)
14 (33.1 mm)	996 (37.7 mm)
15 (35.2 mm)	998 (40.0 mm)

## 53 (36.1 mm) 80 (32.8 mm) 81 (33.2 mm) 82 (35.3 mm) 83 (37.6 mm)

**Candulor TCR Moulds** 

**Candulor NFC+ Moulds**